

Size and colour chart

Glove

SIZES

Glove	Size	a – MCP circumference	b – Wrist circumference
RF001	X SML	16cm – 17.5cm	13.5cm – 15.0cm
RF002	SML	17.5cm – 19.5cm	15.0cm – 16.5cm
RF003	MED	19.5cm – 21.5cm	16.5cm – 18.0cm
RF004	LGE	21.5cm – 23.5cm	18.0cm – 19.5cm
RF005	X LGE	23.5cm – 26.0cm	19.5cm – 21.5cm
RF006	XX LGE	26.0cm – 28.5cm	21.5cm – 23.5cm

COLOURS

Light



Dark




Black



MEASUREMENTS

Measuring Tips:

- All circumferences to be measured to the nearest centimetre.
- Tape should be snug NOT tight
- Position hand as flat as possible



Please note that garments cannot be altered or repaired by Second Skin.



Second Skin Pty Ltd
15/386 Scarborough Beach Road
Osborne Park WA Australia 6017
Email: perth@secondskin.com.au
Phone: +61 8 9201 9455

EC REP

Second Skin
10/32-34 Gordon House Road
Gospel Oak, LONDON NW5 1LP
Email: london@secondskin.com.au
Phone: +44 203 627 9937



Made in Australia

Size and colour chart

Arm Sleeve



By Second Skin

READY TO WEAR GARMENTS

SIZES

Arm Sleeve	Size	a – Upper circumference	b – Elbow circumference	c – Wrist circumference
RF007	X SML	25.5cm – 28.5cm	22.5cm – 24.5cm	13.5cm – 15.0cm
RF008	SML	28.5cm – 31.5cm	24.5cm – 27.0cm	15.0cm – 16.5cm
RF009	MED	31.5cm – 34.5cm	27.0cm – 29.5cm	16.5cm – 18.0cm
RF010	LGE	34.5cm – 37.5cm	29.5cm – 32.5cm	18.0cm – 19.5cm
RF011	X LGE	37.5cm – 41.5cm	32.5cm – 35.5cm	19.5cm – 21.5cm
RF012	XX LGE	41.5cm – 45.5cm	35.5cm – 38.5cm	21.5cm – 23.5cm

COLOURS

Light



Dark




Black



MEASUREMENTS

Measuring Tips:

- All circumferences to be measured to the nearest centimetre.
- Tape should be snug NOT tight
- Arm should be placed at rest by side of the body



Please note that garments cannot be altered or repaired by Second Skin.



Second Skin Pty Ltd
15/386 Scarborough Beach Road
Osborne Park WA Australia 6017
Email: perth@secondskin.com.au
Phone: +61 8 9201 9455

EC REP

Second Skin
10/32-34 Gordon House Road
Gospel Oak, LONDON NW5 1LP
Email: london@secondskin.com.au
Phone: +44 203 627 9937



Made in Australia

Size and colour chart

Knee High Sock Closed Toe

SIZES

Knee High Sock Closed Toe	Size	a – Top circumference	b – Ankle circumference	c – Foot length
RF013	X SML	31.5cm – 33.5cm	28.5cm – 29.5cm	22cm – 23cm
RF014	SML	33.5cm – 35.5cm	29.5cm – 30.5cm	23cm – 24cm
RF015	MED	35.5cm – 37.5cm	30.5cm – 31.5cm	24cm – 25cm
RF016	LGE	37.5cm – 39.5cm	31.5cm – 32.5cm	25cm – 26cm
RF017	X LGE	39.5cm – 41.5cm	32.5cm – 33.5cm	26cm – 27cm
RF018	XX LGE	41.5cm – 43.5cm	33.5cm – 34.5cm	27cm – 28cm

COLOURS

Light



Dark




Black



MEASUREMENTS

Measuring Tips:

- All circumferences and lengths to be measured to the nearest centimetre.
- Tape should be snug NOT tight



Please note that garments cannot be altered or repaired by Second Skin.



Second Skin Pty Ltd
15/386 Scarborough Beach Road
Osborne Park WA Australia 6017
Email: perth@secondskin.com.au
Phone: +61 8 9201 9455

EC REP

Second Skin
10/32-34 Gordon House Road
Gospel Oak, LONDON NW5 1LP
Email: london@secondskin.com.au
Phone: +44 203 627 9937



Made in Australia

Size and colour chart


Knee High Sock Open Toe

SIZES

Knee High Sock Open Toe	Size	a – Top circumference	b – Ankle circumference	c – Foot length
RF019	X SML	31.5cm – 33.5cm	28.5cm – 29.5cm	22cm – 23cm
RF020	SML	33.5cm – 35.5cm	29.5cm – 30.5cm	23cm – 24cm
RF021	MED	35.5cm – 37.5cm	30.5cm – 31.5cm	24cm – 25cm
RF022	LGE	37.5cm – 39.5cm	31.5cm – 32.5cm	25cm – 26cm
RF023	X LGE	39.5cm – 41.5cm	32.5cm – 33.5cm	26cm – 27cm
RF024	XX LGE	41.5cm – 43.5cm	33.5cm – 34.5cm	27cm – 28cm

COLOURS

Light



Dark




Black



MEASUREMENTS

Measuring Tips:

- All circumferences and lengths to be measured to the nearest centimetre.
- Tape should be snug NOT tight



Please note that garments cannot be altered or repaired by Second Skin.



Second Skin Pty Ltd
15/386 Scarborough Beach Road
Osborne Park WA Australia 6017
Email: perth@secondskin.com.au
Phone: +61 8 9201 9455

EC REP

Second Skin
10/32-34 Gordon House Road
Gospel Oak, LONDON NW5 1LP
Email: london@secondskin.com.au
Phone: +44 203 627 9937



Made in Australia

Size and colour chart

Ankle High Sock Closed Toe

SIZES

Ankle High Sock Closed Toe	Size	a – Height from floor	b – Top circumference	c – Ankle circumference	d – Foot length
RF025	X SML	11.0cm	11.0cm – 20.5cm	28.5cm – 29.5cm	22cm – 23cm
RF026	SML	11.5cm	20.5cm – 22.0cm	29.5cm – 30.5cm	23cm – 24cm
RF027	MED	12.0cm	22.0cm – 23.5cm	30.5cm – 31.5cm	24cm – 25cm
RF028	LGE	12.5cm	23.5cm – 24.5cm	31.5cm – 32.5cm	25cm – 26cm
RF029	X LGE	13.0cm	24.5cm – 25.5cm	32.5cm – 33.5cm	26cm – 27cm
RF030	XX LGE	13.5cm	25.5cm – 27.0cm	33.5cm – 34.5cm	27cm – 28cm

COLOURS

Light



Dark




Black



MEASUREMENTS

Measuring Tips:

- All circumferences and lengths to be measured to the nearest centimetre.
- Tape should be snug NOT tight
- Take 'a' measurement to determine location for 'b' measurement



Please note that garments cannot be altered or repaired by Second Skin.



Second Skin Pty Ltd
15/386 Scarborough Beach Road
Osborne Park WA Australia 6017
Email: perth@secondskin.com.au
Phone: +61 8 9201 9455

EC REP

Second Skin
10/32-34 Gordon House Road
Gospel Oak, LONDON NW5 1LP
Email: london@secondskin.com.au
Phone: +44 203 627 9937



Made in Australia

Size and colour chart

Ankle High Sock Open Toe

SIZES

Ankle High Sock Open Toe	Size	a – Height from floor	b – Top circumference	c – Ankle circumference	d – Foot length
RF031	X SML	11.0cm	19.5cm – 20.5cm	28.5cm – 29.5cm	22cm – 23cm
RF032	SML	11.5cm	20.5cm – 22.0cm	29.5cm – 30.5cm	23cm – 24cm
RF033	MED	12.0cm	22.0cm – 23.5cm	30.5cm – 31.5cm	24cm – 25cm
RF034	LGE	12.5cm	23.5cm – 24.5cm	31.5cm – 32.5cm	25cm – 26cm
RF035	X LGE	13.0cm	24.5cm – 25.5cm	32.5cm – 33.5cm	26cm – 27cm
RF036	XX LGE	13.5cm	25.5cm – 27.0cm	33.5cm – 34.5cm	27cm – 28cm

COLOURS

Light



Dark




Black



MEASUREMENTS

Measuring Tips:

- All circumferences and lengths to be measured to the nearest centimetre.
- Tape should be snug NOT tight
- Take 'a' measurement to determine location for 'b' measurement



Please note that garments cannot be altered or repaired by Second Skin.



Second Skin Pty Ltd
15/386 Scarborough Beach Road
Osborne Park WA Australia 6017
Email: perth@secondskin.com.au
Phone: +61 8 9201 9455

EC REP

Second Skin
10/32-34 Gordon House Road
Gospel Oak, LONDON NW5 1LP
Email: london@secondskin.com.au
Phone: +44 203 627 9937



Made in Australia