

# Size and colour chart

## Glove

### SIZES

Glove	Size	a – MCP circumference	b – Wrist circumference
RF001	X SML	16cm – 17.5cm	13.5cm – 15.0cm
RF002	SML	17.5cm – 19.5cm	15.0cm – 16.5cm
RF003	MED	19.5cm – 21.5cm	16.5cm – 18.0cm
RF004	LGE	21.5cm – 23.5cm	18.0cm – 19.5cm
RF005	X LGE	23.5cm – 26.0cm	19.5cm – 21.5cm
RF006	XX LGE	26.0cm – 28.5cm	21.5cm – 23.5cm

### COLOURS

Light



Dark




Black



### MEASUREMENTS

Measuring Tips:

- All circumferences to be measured to the nearest centimetre.
- Tape should be snug NOT tight
- Position hand as flat as possible



Please note that garments cannot be altered or repaired by Second Skin.



Second Skin Pty Ltd  
40 O'Malley Street  
Osborne Park WA Australia 6017  
Email: perth@secondskin.com.au  
Phone: +61 8 9201 9455

EC REP

Second Skin  
10/32-34 Gordon House Road  
Gospel Oak, LONDON NW5 1LP  
Email: london@secondskin.com.au  
Phone: +44 203 627 9937



Made in Australia

# Size and colour chart

## Arm Sleeve

### SIZES

Arm Sleeve	Size	a – Upper circumference	b – Elbow circumference	c – Wrist circumference
RF007	X SML	25.5cm – 28.5cm	22.5cm – 24.5cm	13.5cm – 15.0cm
RF008	SML	28.5cm – 31.5cm	24.5cm – 27.0cm	15.0cm – 16.5cm
RF009	MED	31.5cm – 34.5cm	27.0cm – 29.5cm	16.5cm – 18.0cm
RF010	LGE	34.5cm – 37.5cm	29.5cm – 32.5cm	18.0cm – 19.5cm
RF011	X LGE	37.5cm – 41.5cm	32.5cm – 35.5cm	19.5cm – 21.5cm
RF012	XX LGE	41.5cm – 45.5cm	35.5cm – 38.5cm	21.5cm – 23.5cm

### COLOURS

Light



Dark




Black



### MEASUREMENTS

Measuring Tips:

- All circumferences to be measured to the nearest centimetre.
- Tape should be snug NOT tight
- Arm should be placed at rest by side of the body



Please note that garments cannot be altered or repaired by Second Skin.



Second Skin Pty Ltd  
40 O'Malley Street  
Osborne Park WA Australia 6017  
Email: perth@secondskin.com.au  
Phone: +61 8 9201 9455

EC REP

Second Skin  
10/32-34 Gordon House Road  
Gospel Oak, LONDON NW5 1LP  
Email: london@secondskin.com.au  
Phone: +44 203 627 9937



Made in Australia

# Size and colour chart

## Knee High Sock Closed Toe

### SIZES

Knee High Sock Closed Toe	Size	a – Top circumference	b – Ankle circumference	c – Foot length
RF013	X SML	31.5cm – 33.5cm	28.5cm – 29.5cm	22cm – 23cm
RF014	SML	33.5cm – 35.5cm	29.5cm – 30.5cm	23cm – 24cm
RF015	MED	35.5cm – 37.5cm	30.5cm – 31.5cm	24cm – 25cm
RF016	LGE	37.5cm – 39.5cm	31.5cm – 32.5cm	25cm – 26cm
RF017	X LGE	39.5cm – 41.5cm	32.5cm – 33.5cm	26cm – 27cm
RF018	XX LGE	41.5cm – 43.5cm	33.5cm – 34.5cm	27cm – 28cm

### COLOURS

Light



Dark




Black



### MEASUREMENTS

Measuring Tips:

- All circumferences and lengths to be measured to the nearest centimetre.
- Tape should be snug NOT tight



Please note that garments cannot be altered or repaired by Second Skin.



Second Skin Pty Ltd  
40 O'Malley Street  
Osborne Park WA Australia 6017  
Email: perth@secondskin.com.au  
Phone: +61 8 9201 9455

EC REP

Second Skin  
10/32-34 Gordon House Road  
Gospel Oak, LONDON NW5 1LP  
Email: london@secondskin.com.au  
Phone: +44 203 627 9937



Made in Australia

# Size and colour chart

## Knee High Sock Open Toe

### SIZES

Knee High Sock Open Toe	Size	a – Top circumference	b – Ankle circumference	c – Foot length
RF019	X SML	31.5cm – 33.5cm	28.5cm – 29.5cm	22cm – 23cm
RF020	SML	33.5cm – 35.5cm	29.5cm – 30.5cm	23cm – 24cm
RF021	MED	35.5cm – 37.5cm	30.5cm – 31.5cm	24cm – 25cm
RF022	LGE	37.5cm – 39.5cm	31.5cm – 32.5cm	25cm – 26cm
RF023	X LGE	39.5cm – 41.5cm	32.5cm – 33.5cm	26cm – 27cm
RF024	XX LGE	41.5cm – 43.5cm	33.5cm – 34.5cm	27cm – 28cm

### COLOURS

Light



Dark




Black



### MEASUREMENTS

Measuring Tips:

- All circumferences and lengths to be measured to the nearest centimetre.
- Tape should be snug NOT tight



Please note that garments cannot be altered or repaired by Second Skin.



Second Skin Pty Ltd  
40 O'Malley Street  
Osborne Park WA Australia 6017  
Email: perth@secondskin.com.au  
Phone: +61 8 9201 9455

EC REP

Second Skin  
10/32-34 Gordon House Road  
Gospel Oak, LONDON NW5 1LP  
Email: london@secondskin.com.au  
Phone: +44 203 627 9937



Made in Australia

# Size and colour chart

## Ankle High Sock Closed Toe

### SIZES

Ankle High Sock Closed Toe	Size	a – Height from floor	b – Top circumference	c – Ankle circumference	d – Foot length
RF025	X SML	11.0cm	11.0cm – 20.5cm	28.5cm – 29.5cm	22cm – 23cm
RF026	SML	11.5cm	20.5cm – 22.0cm	29.5cm – 30.5cm	23cm – 24cm
RF027	MED	12.0cm	22.0cm – 23.5cm	30.5cm – 31.5cm	24cm – 25cm
RF028	LGE	12.5cm	23.5cm – 24.5cm	31.5cm – 32.5cm	25cm – 26cm
RF029	X LGE	13.0cm	24.5cm – 25.5cm	32.5cm – 33.5cm	26cm – 27cm
RF030	XX LGE	13.5cm	25.5cm – 27.0cm	33.5cm – 34.5cm	27cm – 28cm

### COLOURS

Light



Dark




Black



### MEASUREMENTS

Measuring Tips:

- All circumferences and lengths to be measured to the nearest centimetre.
- Tape should be snug NOT tight
- Take 'a' measurement to determine location for 'b' measurement



Please note that garments cannot be altered or repaired by Second Skin.



Second Skin Pty Ltd  
40 O'Malley Street  
Osborne Park WA Australia 6017  
Email: perth@secondskin.com.au  
Phone: +61 8 9201 9455

EC REP

Second Skin  
10/32-34 Gordon House Road  
Gospel Oak, LONDON NW5 1LP  
Email: london@secondskin.com.au  
Phone: +44 203 627 9937



Made in Australia

# Size and colour chart

## Ankle High Sock Open Toe

### SIZES

Ankle High Sock Open Toe	Size	a – Height from floor	b – Top circumference	c – Ankle circumference	d – Foot length
RF031	X SML	11.0cm	19.5cm – 20.5cm	28.5cm – 29.5cm	22cm – 23cm
RF032	SML	11.5cm	20.5cm – 22.0cm	29.5cm – 30.5cm	23cm – 24cm
RF033	MED	12.0cm	22.0cm – 23.5cm	30.5cm – 31.5cm	24cm – 25cm
RF034	LGE	12.5cm	23.5cm – 24.5cm	31.5cm – 32.5cm	25cm – 26cm
RF035	X LGE	13.0cm	24.5cm – 25.5cm	32.5cm – 33.5cm	26cm – 27cm
RF036	XX LGE	13.5cm	25.5cm – 27.0cm	33.5cm – 34.5cm	27cm – 28cm

### COLOURS

Light



Dark




Black



### MEASUREMENTS

Measuring Tips:

- All circumferences and lengths to be measured to the nearest centimetre.
- Tape should be snug NOT tight
- Take 'a' measurement to determine location for 'b' measurement



Please note that garments cannot be altered or repaired by Second Skin.



Second Skin Pty Ltd  
40 O'Malley Street  
Osborne Park WA Australia 6017  
Email: perth@secondskin.com.au  
Phone: +61 8 9201 9455

EC REP

Second Skin  
10/32-34 Gordon House Road  
Gospel Oak, LONDON NW5 1LP  
Email: london@secondskin.com.au  
Phone: +44 203 627 9937



Made in Australia